

SkyScan Micro-CT attachment for SEM

Recommended computer configurations

Desktop:

Dell Precision Workstation 490 or T5400:

Two Dual Core Intel® Xeon® 5160 (3.0GHz,1333Mhz,4MB),
Windows® XP Professional, SP2 (NTFS) or Windows Vista 32bit,
4 GB DDR2 667 MHz Quad channel RAM,
128 MB or more, ATI or NVIDIA Graphical card,
300GB or more HDD,
DVD+/-RW Drive,
19"or more LCD UltraSharp Monitor,
Optical Mouse with wheel,
Internal speaker or AS501 SoundBar.

Laptop:

Sony VAIO VGN-FZ21Z:

Intel® Core™ 2 Duo Processor T7500 (2.2GHz,800Mhz,4MB),
Windows Vista 32bit,
2 GB DDR2 667 MHz SDRAM,
256MB-1023MB NVIDIA GeForce 8600M Graphics,
300GB HDD,
CD/DVD/ Blu-ray +/-RW Drive,
15.4" LCD X-black LCD with double lamp technology,
Internal speakers and touchpad.
Additional optical mouse with wheel recommended.